

Operating Manual

MASTERPULS MP100



Inhalt

1	General Information	5
1.1	Introduction	5
1.1.1	Indications	6
1.1.2	Contraindications	7
1.1.3	Side effects	7
1.2	Key	8
1.3	Description of controls and functional elements	9
1.3.1	MASTERPULS MP 100	9
1.3.2	Applicator	11
1.3.3	Optional ENERGY compressor for MASTERPULS MP100 (Optional)	13
2	Installation	14
2.1	Unpacking	14
2.2	Scope of supply	14
2.3	Installation	15
2.3.1	Applicator holder installation	15
2.3.2	Mains voltage selection	16
2.3.3	Connecting the compressed air supply	17
2.3.4	Connecting power supply cables between systems and socket	19
2.3.5	Applicator connection	19
2.3.6	Optional foot switch connection (Foot switch KARL STORZ 20010230)	20
2.3.7	Optional potential equalization	20
3	Operation	21
3.1	General warnings and safety information	21
3.2	Start up	23
3.3	Functional checks	25
3.4	Standard settings	25
3.5	Treatment	26
4	Cleaning, maintenance, overhaul	27
4.1.1	Control unit	27
4.1.2	Applicator	27
4.1.3	Optional foot switch (KARL STORZ 20010230)	28
4.1.4	Optional ENERGY compressor	28
4.1.2.1	Emptying the condensate collector	28
4.1.2.2	Replacing the filter element	29
4.1.4.2	ENERGY compressor mains fuse replacement	33
4.2	Maintenance	34
4.3	Disposal	34
4.4	Repair	34
4.4	Applicator sealing ring replacement	35

TABLE OF CONTENTS

4.5	Replacement of shock absorber in applicator	36
4.6	Overhaul	36
4.6.1	Content of the overhaul kit	37
4.6.2	Overhaul instructions	37
5	Trouble-shooting	41
5.1	Trouble-shooting	41
6	Accessories and spare parts	42
6.1.1	MASTERPULS MP100	42
6.1.2	Optional ENERGY compressor	43
6.1.3	Accessories	43
6.1.4	Documentation	43
7	Technical specifications	44
7.1	MASTERPULS MP100	44
7.2	Compressor	45
7.2.1	Energy compressor	45
7.3	Conformity with standards	46
7.2.2	Sil.Air 50 TDC compressor	46
7.4	Conformity with directives	47
8	Warranty and Service	47
8.1	Warranty	47
8.2	Service	47

1 General Information

1.1 Introduction

This manual contains warnings, safety instructions and specific operating instructions in accordance with liability regulations.

CAUTION

Complete or partial failure to observe the instructions, information or procedures preceded by the term „CAUTION“ may cause injury or fatal accidents.



ATTENTION

Complete or partial failure to observe the instructions, information or procedures preceded by the term „ATTENTION“ may cause equipment damage.



NOTE

Additional information concerning specific features or operating instructions is preceded by the term „NOTE“.

CAUTION

Before you start using the MASTERPULS MP100 for the first time, please make sure you have read and understood all information provided in this operating manual.

Please refer to the separate operating manual prior to any use or overhaul of the Sil.Air 50 TDC compressor.

Familiarity with the information and instructions contained in this manual is an essential requirement to ensure efficient and optimal use of the system, to avoid dangers to persons and to the equipment and to obtain good treatment results.

Thorough knowledge of the information included in this manual will also enable you to react promptly and effectively in case of malfunctions and failures.



The MASTERPULS MP100 is a compressed air operated ballistic shock wave generator featuring high-precision ballistic components in its applicator for shock wave generation. The motion and weight of the projectile accelerated by compressed air produce kinetic energy that is converted into sound energy when the projectile strikes an unmoved surface (shock transmitter). This acoustic pulse is transmitted to the target tissue by means of an acoustic impedance adapter (shock wave coupling cushion) or gel.

NOTE

Medical devices operating on the basis of the above principle are generally referred to as extracorporeal shock wave systems in modern medical literature.

1.1.1 Indications

The MASTERPULS MP100 is designed for extracorporeal shock wave therapy. Indications include:

- calcaneal spur / plantar fascitis
- shoulder pain with or without calcifications
- achillodynia
- trochanteric bursitis / proximal iliotibial band friction syndrome
- radial/ulnar humeral epicondylitis
- patellar tendinopathy
- tibial edge syndrome
- insertion tendonitis in general
- acupuncture
Qualified training in acupuncture and acupuncture shock wave therapy (AkuST) is required for therapeutic application of the MASTERPULS MP100 in the field of acupuncture.
- muscle trigger points
A sound knowledge of trigger point therapy and trigger point shock wave therapy (TrST) is required for therapeutic application of the MASTERPULS MP100 in the field of trigger point shock wave therapy.

1.1.2 Contraindications

Treatment with the STORZ MEDICAL MASTERPULS MP100 is not permitted in the following cases:

- coagulation disorders (haemophilia)
- use of anticoagulants, especially Marcumar
- thrombosis
- tumour diseases, carcinoma patients
- pregnancy
- polyneuropathy in case of diabetes mellitus
- acute inflammations / pus focus in the target area
- children in growth
- cortisone therapy up to 6 weeks before first treatment

CAUTION

Shock waves must not be applied to target areas located above air filled tissue (lungs), nor to any regions near large nerves, vessels, the spinal column or head.
















1.1.3 Side effects

Treatment with the MASTERPULS MP100 may cause the following side effects:

- swelling, reddening, haematomas
- petechiae
- pain
- skin lesions after previous cortisone therapy

These side effects generally abate after 5 to 10 days.

1.2 Key

	Operating manual to be observed!
	Applicator connector
	Foot switch connector (optional)
	Potential equalization
	AC
I/O	ON/OFF switch
	Application pressure in bar (pressure gauge reading)
	Analog setting of application pressure (control knob)
M	Mode button for mode selection: single shock / continuous shock
	Single shock
	Continuous shocks at 5 Hz frequency Continuous shocks at 15 Hz frequency / blinking
	Continuous shocks at 10 Hz frequency
	Applied shocks in total (digital display)
	Fuse
	Type B applied part
Air: 5-7bar	Compressed air connector

1.3 Description of controls and functional elements

1.3.1 MASTERPULS MP 100




Abb. 1 - 1 Front view of MASTERPULS MP100

- 1 Pressure gauge indicating selected energy level (application pressure)
- 2 Selected shock frequency symbol: single shock
- 3 Selected shock frequency symbol: 5 Hz, 10 Hz, 15 Hz
- 4 Treatment shock counter (resettable)
- 5 Shock counter reset button
- 6 Foot switch connector
- 7 Applicator connector
- 8 Mode button, shock frequency selector button
- 9 Energy level control knob (application pressure)
- 10 Mains switch

CHAPTER 1 - GENERAL INFORMATION




Fig. 1 – 2 Rear view of MASTERPULS MP100

- 1 Operation shock counter (not resettable)
- 2 Mains connector
- 3 Fuse holder
- 4 Potential equalization connector
- 5 Compressed air connector (input)

1.3.2 Applicator




Fig. 1 – 3 Applicator with connecting cable

- 1 Shock transmitter
- 2 Shock transmitter screw cap
- 3 Trigger button
- 4 Locking button
- 5 Shock absorber
- 6 Applicator connecting cable




Fig. 1 – 4 Shock transmitter with sealing rings

CHAPTER 1 - GENERAL INFORMATION

Depending on the therapy to be performed, the applicator can be equipped with one of the following five shock transmitters:

- ESWT shock transmitter (Fig. 1 – 5/1) for ESWT therapy: Ø 15 mm
- TrST shock transmitter (Fig. 1 – 5/2) for trigger point therapy: Ø 10 mm
- AkuST shock transmitter (Fig. 1 – 5/3) for therapeutic shock wave application in acupuncture: Ø 6 mm (Fig. 1 - 5/3) Ø 6mm




Fig. 1 – 5 ESWT, TrST and AkuST shock transmitters

- Optional: Focus shock transmitter
- Optional: D-Actor shock transmitter

NOTE

Always use an acoustic impedance adapter (shock wave coupling cushion) (Fig. 1 – 6/1) in order to ensure perfect hygiene during patient treatment.



Fig. 1 – 6 Shock transmitter with acoustic impedance adapter

1.3.3 Optional ENERGY compressor for MASTERPULS MP100 (Optional)

The ENERGY compressor provides the compressed air supply for the MASTERPULS MP100 and can be ordered as an option.




Fig. 1 – 7 Front view of ENERGY compressor

1 Mains switch




Fig. 1 – 8 Rear view of ENERGY compressor

- 1 Compressed air connector (output)
- 2 Fuse holder
- 3 Mains connector
- 4 Connector for mains power link cable (mains output for MASTERPULS MP100 control system)
- 5 Potential equalization connector
- 6 Condensate collector

2 Installation

2.1 Unpacking

- Remove the equipment and accessories from the packaging container. Proceed with extreme caution.
- Check that all items are included in the packaging container and that they are not damaged.
- Contact your supplier or the manufacturer immediately if any items are missing or damaged. Retain the original packaging, if possible. It may prove useful for any later equipment transport.

2.2 Scope of supply

The standard scope of supply of the STORZ MEDICAL MASTERPULS MP100 includes the following items:

- MASTERPULS MP100 (control system)
- applicator with accessories (2 shock transmitters, complete applicator holder, acoustic impedance adapter)
- mains cable
- mains power link cable
- compressed air tube
- overhaul kit for applicator overhaul
- 850 ml gel
- user manual (operating manual, system logbook and training records)



Please refer to Section 6 for information on ACCESSORIES AND SPARE PARTS.

2.3 Installation

2.3.1 Applicator holder installation

The applicator can be fitted to the right or left side of the system as desired by the system user. Owing to the specific shape of the applicator holder and the location of the bores provided on the side panels of the system, the applicator holder can be mounted either parallel to the system front (Fig. 2 – 1/1) or in an inclined position.

- If you wish to mount the applicator holder parallel to the system front, fix the holder by tightening both screws (Fig. 2 – 1/).



Fig. 2 – 1 Mounting position parallel to system front

- If you wish to mount the applicator holder in such a way that it is inclined towards the front, then use only one screw (Fig. 2 – 2) to fix the holder, making sure that the bevelled edge of the holder ends at the front strip of the side panel.



Fig. 2 - 2 Inclined mounting position

- Regardless of the installation position you choose, make sure that the smaller holder ring points downwards.

2.3.2 Mains voltage selection



ATTENTION

Check that the mains voltage set on the MASTERPULS MP100 matches the local mains voltage to avoid system damage.

- Remove the mains fuse holder (Fig. 2 – 3/1) from the rear panel of the MASTERPULS MP100 and take off the voltage selector (Fig. 2 – 3/3)




Fig. 2 – 3 Mains fuse holder

NOTE

The MASTERPULS MP100 is set to a local mains voltage of 230 VAC when delivered.

- Make sure to use only fuses that meet the requirements specified in Table 2 – 1.
- Replace the fuses, if necessary.

Mains voltage	Fuse rating
230 / 240 VAC	2 x T 0,2 A L
100 / 120 VAC	2 x T 0,5 A L


Table 2 – 1 Fuse rating

- Insert the voltage selector into the mains fuse holder in such a way that the mains voltage visible through the glass window on the outside of the fuse holder (Fig. 2 – 3/2) matches the local mains voltage.

ATTENTION

Check that the fuses installed are rated for the selected voltage!

The mains voltage of the optional compressor cannot be changed. The compressor is available in three different versions (110 V, 115 V, 220 - 240 V).



2.3.3 Connecting the compressed air supply

ATTENTION

The compressed air input pressure must be 5.0 - 7.0 bar. If compressors other than the STORZ MEDICAL ENERGY compressor are used, this pressure must be checked on the compressed air supply system.

If the pressure exceeds 7.0 bar the system's overpressure protection will be tripped.

The compressed air must be dry, oil-free and filtered to ≤ 5 micron.

Never use CO2 gas for compressed air supply.





ATTENTION

If compressed air is supplied through a wall outlet, use a pressure reducer (max. 7 bar), if necessary.

- Place the optional ENERGY compressor onto a solid and stable surface. The MASTERPULS MP100 and the ENERGY compressor may be positioned on top or beside each other.
- Connect the compressed air tube to the compressed air connectors provided on the compressor (Fig. 2 – 5/1) and MASTERPULS MP100 (Fig. 2 – 4/5).




Fig. 2 – 4 Connectors on the MASTERPULS MP100




Fig. 2 - 5 Connectors on the ENERGY compressor

2.3.4 Connecting power supply cables between systems and socket

- Connect the mains power link cable to the mains output connector (Fig. 2 – 5/3) provided on the ENERGY compressor.
- Connect the other end of the mains power link cable to the mains connector of the MASTERPULS MP100 (Fig. 2 – 4/2).
- Connect one end of the mains cable to the mains connector (Fig. 2 – 5/2) of the compressor and the other end to the socket.
- If the MASTERPULS MP100 is used without the optional ENERGY compressor, the mains cable is directly connected to the MASTERPULS MP100 and to the socket.

2.3.5 Applicator connection

- Connect the applicator cable to the applicator connector (Fig. 2 – 6/2) provided on the MASTERPULS MP100. Make sure that the red spots on the cable connector match the red spots on the applicator connector.




Fig. 2 – 6 Applicator connector

2.3.6 Optional foot switch connection (Foot switch KARL STORZ 20010230)

- Connect the foot switch cable (special accessories) to the foot switch connector (Fig. 2 – 6/1) provided on the MASTERPULS MP100.

NOTE

The optional foot switch (KARL STORZ 20010230) complies with IPX8 rating according to IEC529 (degree of protection against water).

2.3.7 Optional potential equalization

The MASTERPULS MP100 and the ENERGY compressor both feature a potential equalization connector (see Fig. 2 – 4/4 and Fig. 2 – 5/4). Where necessary, connections for potential equalization must be made by suitably qualified personnel.

3 Operation

3.1 General warnings and safety information

CAUTION

The MASTERPULS MP100 must only be used by suitably qualified and trained medical personnel.

The MASTERPULS MP100 must only be used for applications approved by STORZ MEDICAL AG!

To avoid safety hazards, use of the system for applications other than those specified in Section 1.1.1 INDICATIONS is not allowed!

Unplug the mains cable from the system and compressor before carrying out cleaning and maintenance work!

Do not use the MASTERPULS or the ENERGY compressor in potentially explosive environments, i.e. in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

The optional KARL STORZ foot switch must not be used in potentially explosive atmospheres according to classification AP as per IEC 60601.

Unplug the applicator cable from the control system before carrying out cleaning and maintenance work. Do not reconnect the applicator until it has been completely reassembled!

Do not try to open the system! Risk of electric shocks!

Risk of transmission of microorganisms! Disinfect the applicator after each use! Refer to chapter 4 CLEANING for details.



→ 1.1.1

→ 4



ATTENTION

Check that the installation surfaces have sufficient carrying capacity to avoid equipment damage!

The system must only be connected to properly earthed and correctly installed shockproof sockets!

Check that the system is in perfect working order before each use.

Never cover the system when in use!

Make absolutely sure that no liquid can seep into the system housing or applicator.

**ATTENTION**

Equipment damage caused by failure to use the system in accordance with its intended use or by operator errors is not covered by the warranty.

The MASTERPULS MP100 must only be used with accessories that have been approved by the system manufacturer. To prevent safety hazards, unauthorized system modifications are not allowed and will void the CE mark approval and warranty.

Disposal of the system and its components must be carried out in accordance with national waste disposal regulations.

NOTE

The MASTERPULS MP100 as well as the ENERGY compressor meets the requirements of the applicable electromagnetic compatibility (EMC) standards EN 60601-1-2.

These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected
- Consult the manufacturer of field service technician for help.

3.2 Start up

- Switch on the compressor followed by the control system.
- Use the control knob (Fig. 3 – 1/5) to select an initial value for the shock application pressure (min. 1.0 bar). The selected value is indicated on the pressure gauge (Fig. 3 – 1/1).




Fig. 3 – 1 Application pressure selection on control system

The maximum application pressure is limited to 4 bar. To ensure correct system operation, a minimum 1.0 bar pressure must be selected.

- Select the shock frequency by pressing the M (mode) button (Fig. 3 – 1/4).



single shock



continuous shocks: 5 shocks / sec.



continuous shocks: 10 shocks / sec.



continuous shocks: 15 shocks / sec., press the mode button for at least 2 sec.

- Press the trigger button provided on the applicator (Fig. 3 – 2/2).
- A single press of the trigger button will cause one shock to be released.
- In order to work in continuous shock mode, press and hold the trigger button or press the locking button (Fig. 3 – 2/1) while at the same time pressing the trigger button.

CHAPTER 3 - OPERATION

- Unlock the locking button if you wish to disable continuous shock mode.




Fig. 3 – 2 Trigger button and locking button

3.3 Functional checks

Perform the following functional checks after the system has been installed:

- Check the control system and applicator for damage.
- Start the MASTERPULS MP100.
- Set the application pressure to 2.0 bar.
- Reset the treatment shock counter on the front of the control system.
- Release individual shocks in single shock mode.
- Release shocks in continuous shock mode (shock frequency 5 Hz , 10 Hz and 15 Hz).
- Check that the triggered shocks are correctly counted on the treatment shock counter on the front of the control system.
- Select the maximum 4 bar application pressure.
- Release individual shocks in single shock mode.
- Release shocks in continuous shock mode (shock frequency 5 Hz , 10 Hz and 15 Hz).
- Release shocks by means of the foot switch, if used.

3.4 Standard settings

- Before each therapy session, press the reset button (Fig. 3 – 1/3) to reset the treatment shock counter (Fig. 3 – 1/2) provided on the control system.
- The number of shocks to be applied for trigger point shock wave therapy (TrST) and acupuncture shock wave therapy (AkuST) differs from ESWT treatment. Such therapies must only be performed by suitably qualified and trained personnel.
- A total of about 2000 shocks must generally be applied per therapy session. Please refer to the MASTERPULS MP100 application brochure for details.
- Start the therapy at a pressure of 2 bar and a frequency of 5 Hz.

3.5 Treatment

- Perform treatment as described in the MASTERPULS MP100 application brochure.
- Apply a sufficient amount of coupling gel to the patient's skin in the coupling area and to the shock transmitter.
- When using the acoustic impedance adapter, generously apply coupling gel to both sides of the adapter. Bear in mind that the acoustic impedance adapter cannot be used if treatment is performed with the TrST shock transmitter.
- Fit the acoustic impedance adapter into the mounting groove (Fig. 3 – 3/1) provided on the shock transmitter screw cap.




Fig. 3 – 3 Mounting groove for acoustic impedance adapter



➔ 1.1.1

CAUTION

The MASTERPULS MP100 must only be used for applications approved by STORZ MEDICAL AG!

To avoid safety hazards, use of the system for applications other than those specified in Section 1.1.1 (INDICATIONS) is not allowed!

4 Cleaning, maintenance, overhaul

4.1 Cleaning

Regular cleaning of the system ensures perfect hygiene and operation of the MASTERPULS MP100.

CAUTION

Prior to performing any cleaning unplug the mains cable.



4.1.1 Control unit

- Clean the exterior of the unit with a moist cloth. Use soap water or a mild cleaning agent.

ATTENTION

Under all circumstances the unit tubes must be protected from intruding fluids of all kinds.



4.1.2 Applicator

The applicator must be thoroughly cleaned and disinfected after each therapy session.

- Remove the acoustic impedance adapter, if necessary.
- Clean the applicator with soap water to remove the ultrasound coupling gel.
- We advise to open the shock transmitter screw cap and to clean the shock transmitter as well as the screw cap separately.
- Disinfect the applicator with a surface disinfectant.
- The acoustic impedance adapter is not suitable for reuse.

4.1.3 Optional foot switch (KARL STORZ 20010230)

- Clean the foot switch with soap water or a mild cleaning agent.

NOTE

The foot switch is classified as IPX8 (degree of protection against water) according to IEC529.

4.1.4 Optional ENERGY compressor

- Clean the exterior of the compressor with a wet cloth. Use soap water or a mild cleaning agent.



ATTENTION

Under all circumstances the unit tubes must be protected from intruding fluids of all kinds.

4.1.2.1 Emptying the condensate collector

The condensate collector on the rear side of the ENERGY compressor should be emptied when three quarters full or earlier.

- Remove the condensate collector (Fig. 4 – 1/1) from the holder on the rear side of the compressor.
- Empty the collector and place it back into the holder.




Fig. 4 – 1 Condensate collector

4.1.2.2 Replacing the filter element

CAUTION

The pump head of the ENERGY compressor becomes very hot during operation. Let the compressor cool down before replacing the filter element.



Contents of the filter element replacement kit (Art.-No. 14464):

- 1 filter element
- 1 allen key
- 1 air input filter

The filter element and the air input filter of the ENERGY compressor should be replaced in case of performance loss (drop of pressure during the release of shocks). Proceed as follows:

- Switch off the compressor by setting the mains switch provided on the front panel to OFF. Unplug the mains cable.
- Vent the compressor system before removing the compressed air tube. To do so, press the rapid action safety coupling (Fig. 4 – 2/1) provided for the compressed air tube towards the rear panel of the compressor until it snaps into place.




Fig. 4 – 2 System venting

- After having vented the compressor, pull the rapid action safety coupling back into its normal position.
- Remove the compressed air tube from the rapid action safety coupling.
- Turn the compressor upside down in order to be able to open the bottom panel.
- Remove the 4 hexagon socket screws (Fig. 4 – 3/1 to Fig. 4 – 3/4) from the bottom panel by means of the supplied allen key.

CHAPTER 4 - CLEANING, MAINTENANCE, OVERHAUL




Fig 4 – 3 Opening the bottom panel

- Remove the condensate discharge tube (Fig. 4 – 4/1) from the L-type screwed tube joint of the metal container. While doing so, press the metal clamping ring (Fig. 4 – 4/2) backwards.




Fig. 4 – 4 Opening the metal container

- Unscrew the metal container (Fig. 4 – 4/4) from its mounting, eventually using a pair of pliers. For this purpose, please seize the nut of the L-type screwed tube joint (Fig. 4 – 4/3).
- The pressure filter metal container (Fig. 4 – 5/1) is equipped with a float valve (Fig. 4 – 5/2) and sealing ring (Fig. 4 – 5/3). Make sure that during removal of the metal container the float valve and sealing ring remain attached to the container and do not fall into the compressor.

CHAPTER 4 - CLEANING, MAINTENANCE, OVERHAUL




Fig. 4 – 5 Metal container with float valve and sealing ring

- Screw off the white filter element (Fig. 4 – 6/1) with the black air current control ring (Fig. 4 - 6/2).



Fig. 4 – 6 Removal of filter elements

- Screw the air current control ring (Fig. 4 – 7/2) off the white filter element. Dispose of the filter element (Fig 4 – 7/1). The air current control ring will be reused.




Fig. 4 - 7 White filter element with current control ring

- Take the new white filter element included in the filter element replacement kit and screw it into the air control ring .
- Install the filter element in the bottom of the pressure filter metal container (Fig. 4 – 6/3) in the compressor.

- If the float valve (Fig. 4 – 8/2) falls out of the new pressure filter metal container, move it back into its original position by means of tweezers so that the sealing pin (Fig. 4 – 8/1) is located precisely in the opening provided inside the metal container.



Fig. 4 - 8 Float valve with sealing pin



ATTENTION

If the sealing pin is bent or if it is positioned inaccurately, the float valve installed in the pressure filter metal container may come loose and may not properly seal off any more!

- Insert the new sealing ring provided for the pressure filter metal container (Fig. 4 – 5/3) into the cavity on the bottom edge of the container. Check that the sealing ring fits well. If necessary, moisten the sealing ring to make sure it cannot fall out.
- Screw the pressure filter metal container to the holder. Tighten it until hand-tight.
- Push the condensate discharge tube into the L-type screwed tube joint.
- Unscrew the used air input filter (Fig. 4 - 4/5), preferable using a pair of pliers, and screw the new air input filter part of the filter element replacement kit on to mounting.
- Fasten the bottom panel by means of the 4 screws
- Fasten the compressed air tube in the rapid action safety coupling.
- Plug in the mains cable.

4.1.4.2 ENERGY compressor mains fuse replacement

The holder of the mains fuse is located at the rear panel of the compressor (see Fig. 1 - 8/2).

- Push the clip of the mains fuse holder (Fig. 4 - 9/1) upwards and take the holder off the housing.




Fig. 4 - 9 Mains fuse holder

- Pull the old fuses out of the fuse holder (Fig. 4 - 10/1).




Fig. 4 - 10 Fuse replacement

- Replace the fuses.
- Push the fuse holder back in to the opening until engaged.

4.2 Maintenance

Preventive maintenance is not necessarily required. However, regular maintenance may help to identify possible defects at an early stage and thus increase the safety and service life of the equipment.

Maintenance services can be ordered from our regional representatives in your area or directly from STORZ MEDICAL AG.

We recommend that functional and safety checks be performed at least once a year. National accident prevention regulations and test and inspection intervals prescribed for medical devices must, of course, be observed.

NOTE

For further details on content and performance of the safety checks please contact your local dealer.

Perform the following tests in order to ensure perfect functioning of the MASTERPULS MP100 and the ENERGY compressor.

1. Leakage current test at chassis < 100 μ A
2. Earth impedance test (incl. applicator housing) < 0,2 Ohm (with mains cable)
3. Power consumption test \leq rated power

For additional information on the overhaul of the Sil.Air 50 TDC compressor please refer to the separate operating manual.

4.3 Disposal

When disposing of the present medical device products no special measures have to be observed. Please proceed in accordance with the national regulations. After expiration of their life time, dispose the MASTERPULS MP100 and the ENERGY compressor as electronical scrap.

Additional information about the disposal of the Sil.Air 50 TDC compressor please find in the separate operating manual.

4.4 Repair

Repair work on defective equipment must only be carried out by personnel suitably authorized by STORZ MEDICAL. Only original STORZ MEDICAL spare parts must be used for this purpose.

4.4 Applicator sealing ring replacement

NOTE

The applicator is equipped with a projectile safety catch device in order to prevent the projectile from exiting in case a shock is released while shock transmitter as well as shock transmitter screw cap are taken off. The safety catch is also activated when the shock transmitter screw cap is not well tightened or the sealing ring between screw cap and shock transmitter is missing, respectively at the shock transmitter endings two sealing rings old/new come together.

Once the projectile gets caught by the safety catch, it can be pushed back into the guide tube using reasonable force.

If this is repeatedly the case or if it is impossible to release the projectile or if the receptory edge is defectuous, the applicator shaft should be replaced .

The applicator sealing rings must be replaced after about 250,000 shocks or earlier if leaks are identified on the applicator. Proceed as follows:

- Switch off MASTERPULS MP100.
- Unplug the applicator cable from the MASTERPULS MP100.
- Screw off the shock transmitter screw cap (Fig. 4 – 9/1).
- Remove both sealing rings (Fig. 4 – 9/2) and dispose of them.

ATTENTION

Make sure that no worn sealing rings are in the screw cap or applicator!

- Clean shock transmitter and shock transmitter screw cap from remaining oil.
- Fit two new sealing rings to the shock transmitter (Fig. 4 – 9/5).
- Place the shock transmitter with the new sealing rings into the shock transmitter screw cap.
- Firmly screw the shock transmitter screw cap to the applicator.






Fig. 4 – 11 Sealing ring replacement

4.5 Replacement of shock absorber in applicator

- Remove the old shock absorber (Fig. 4 – 11/4) with a pointed object (screwdriver).
- If necessary, remove remaining adhesive elements.
- Take the new shock absorber included in the applicator overhaul kit and remove the protective film from the two adhesive strips.
- Push the new shock absorber into the cavity on the applicator handle until it snaps into place in the cavity frame.

4.6 Overhaul

Shock waves are generated mechanically. Due to the effects of friction, the applicator components are continuously exposed to mechanical stress that will cause minor wear.

NOTE

The applicator should be overhauled about every 750,000 shocks. Using the overhaul kit that includes all required wear parts, the applicator can be overhauled easily and quickly by the user of the equipment.

4.6.1 Content of the overhaul kit

Article no.	Product description
13692	ESWT shock transmitter (Ø 15 mm) with sealing rings
13280	Projectile
13271	Guide tube
15558	Sealing rings MP100
13677	Fork wrench (jaw size 22)
13479	Box with protective foam sheets
13272	Shock absorber

The overhaul kit can be ordered from your dealer quoting order number 13460.

4.6.2 Overhaul instructions

NOTE

The applicator is equipped with a projectile safety catch device in order to prevent the projectile from exiting in case a shock is released while shock transmitter as well as shock transmitter screw cap are taken off. The safety catch is also activated when the shock transmitter screw cap is not well tightened or the sealing ring between screw cap and shock transmitter is missing, respectively at the shock transmitter endings two sealing rings old/new come together.

Once the projectile gets caught by the safety catch, it can be pushed back into the guide tube using reasonable force.

If this is repeatedly the case or if it is impossible to release the projectile or if the receptory edge is defectuous, the applicator shaft should be replaced .

- Switch off the MASTERPULS MP100.
- Unplug the applicator cable from the control system.
- Place the applicator onto a dry, clean and dust-free surface.

CHAPTER 4 - CLEANING, MAINTENANCE, OVERHAUL

- Screw off the shock transmitter screw cap (Fig. 4 – 12/1).



Fig. 4 – 12 Applicator overhaul

- Using the supplied fork wrench (Fig. 4 – 13/1), remove the shaft (Fig. 4 – 13/2) from the applicator and pull it out of the applicator handle (Fig. 4 – 13/3).




Fig. 4 – 13 Shaft removal

- Remove the guide tube (Fig. 4 – 14/1) from the shaft.




Fig. 4 – 14 Guide tube removal

CHAPTER 4 - CLEANING, MAINTENANCE, OVERHAUL

- A magnet is provided in the applicator handle (Fig. 4 – 13/3) to hold the projectile. To remove the projectile, hold the applicator handle in such a way that its opening points downwards. Slightly knock the handle against the table surface until the projectile falls out.
- In case the projectile breaks apart due to overloading, projectile fragments may be located inside the guide tube.
- Dispose of the worn guide tube and of all worn shock transmitters and sealing rings.
- Clean the shaft and shock transmitter screw cap with alcohol. These parts will be reused.
- Make sure to place all cleaned components onto a dry, clean and dust-free surface.
- Take the new guide tube, projectile and sealing rings included in the overhaul kit.
- Insert the guide tube (Fig. 4 – 15/3) into the opening in the applicator handle (Fig. 4 – 15/1) up to the limit stop.




Fig. 4 – 15 Guide tube and projectile fitting

- Make sure that the guide tube end inserted into the applicator handle is the one featuring two air slots (Fig. 4 – 16/1).




Fig. 4 – 16 Guiding tube air slots



→ 3.3

- Insert the new projectile (Fig. 4 - 15/3) into the fitted guide tube.
- Screw the shaft into the applicator and tighten it by means of the supplied fork wrench until hand-tight.
- New shock transmitters come with new sealing rings. Thus, only replace the sealing ring on the shaft (Fig. 4 - 12/2).

ATTENTION

Make sure that no worn sealing rings are in the shock transmitter screw cap and on the shaft end!

- Replace the shock absorber.
- Screw the screw cap with the required shock transmitter back onto the shaft.
- Carry out a functional check of the applicator as described in Section 3.3 (FUNCTIONAL CHECKS).

5 Troubleshooting

5.1 Troubleshooting

CAUTION

Unplug the mains cable from the system before you carry out any maintenance work!



Fault	Possible cause	Corrective actions
System does not work Operating mode LED not lit	Power failure	Check the power supply.
	Defective mains fuse	Replace the fuses.
	Defective mains plug	Replace the mains cable.
No compressed air supply	Defective foot switch connecting cable	Check the cable and tube connections and replace them, if necessary.
	Leaks on applicator cable or cable not properly connected	Check the compressor air filter and replace it, if necessary.
	Compressed air tube not connected or not correctly fastened with safety coupling	
	Clogged compressor air filter	
No shock wave power output	No compressed air supply	Check the compressed air supply.
	Blocked or worn projectile	Take applicator apart. Clean the guide tube and projectile.
		Overhaul the applicator.

6 Accessories and spare parts

6.1.1 MASTERPULS MP100

Holder, complete	13701
------------------	-------

Compressed air tube

Compressed air tube for compressor, 1 m long	13463
--	-------

Compressed air tube for compressor, 3 m long	13447
--	-------

Compressed air tube for central compressed air supply	13464
---	-------

Mains cable

Mains cable CEE 7 Europe, 4 m long (for control system or SMAG compressor)	13455
---	-------

IEC coupling, 1 m long (between control system and compressor)	13546
---	-------

Applicator

Ø 15 mm ESWT shock transmitter set	13476
------------------------------------	-------

Ø 10 mm TrST shock transmitter set	13457
------------------------------------	-------

Ø 6 mm AkuST shock transmitter set	13458
------------------------------------	-------

Ø 20 mm D-Actor shock transmitter set	13459
---------------------------------------	-------

Ø 15 mm Focus shock transmitter set	15100
-------------------------------------	-------

Applicator shaft	14615
------------------	-------

Applicator shaft with notch	15236
-----------------------------	-------

Overhaul

Applicator overhaul kit	13460
-------------------------	-------

D- actor overhaul kit	13477
-----------------------	-------

6.1.2 Optional ENERGY compressor

Filter replacement set	14464
Insulating pad	14463

6.1.3 Accessories

Gel bottle set (10 bottles)	13473
MASTERPULS MP100 transport suitcase (for control system or Energy compressor or Sil.Air 50 TDC compressor)	13472
MASTERPULS MP100 transport suitcase (for control system or compressor)	11608
MASTERPULS MP100 Trolley	13471
MASTERPULS MP100 Trolley with compressor covering	15149
Foot switch (KARL STORZ 20010230) incl. connection cable	10103
Acoustic impedance adapter	13470
F- Meter	15323

6.1.4 Documentation

MASTERPULS MP100 user manual German/English/French	13833
---	-------

7 Technical specifications

7.1 MASTERPULS MP100

Operating modes	single shock, continuous shock 5 Hz / 10 Hz / 15 Hz
Energy selection	continuously adjustable
Max. sound pulse amplitude	14,3 MPa
Mains input voltage	100 / 120 / 230 / 240VAC $\sim \pm 10\%$
Mains frequency	50 / 60Hz
Mains fuse	230 / 240VAC : 2 x T 0,2 A L 100 / 120VAC : 2 x T 0,5 A L
Power consumption	40 VA
Compressed air supply	5 - 7 bar
Compressed air output	Max. 4 bar
Ambient temperature during operation	10 – 40°C
Ambient temperature during storage and transport	0 – 60°C
Ambient air pressure	500 hPa - 1060 hPa
Air humidity	5 – 95%, non condensing
Control system weight	5,2 kg
Applicator weight	0,44 kg
Housing dimensions (wxhxd)	305 x 155 x 233 mm
Classification	Class I equipment
Protection against ingress of water	IPX0

This device complies with the standards EN 60601-1, C22.2 Nr. 601.1-M90 and UL std N0 60601-1.

7.2 Compressor

7.2.1 Energy compressor

Art. No.	13465	13466	13476
Mains voltage	220-230 V	115 V	100 V
Mains frequency	50 / 60 Hz	60 Hz	50 / 60 Hz
Mains fuse	T 5 A L (SB)	T 6,3 A L (SB)	T 6,3 A L (SB)
Compressed air output	6 bar		
Power consumption	max. 1225 VA		
Power consumption: compressor without additional consumer	450 VA	500 VA	500 VA
Ambient temperature during operation	10 – 40°C		
Ambient temperature during storage and transport	0 – 60°C		
Ambient air pressure	500 hPa - 1060 hPa		
Air humidity	5 – 95%, non-condensing		
Weight	10,9 kg	10,2 kg	10,2 kg
Housing dimensions (wxhxd)	305 x 155 x 233 mm		
Classification	class I equipment		
Protection against ingress of water	IPX0		

This device complies with the standards EN 60601-1, art. No 13466 complies also with C22.2 Nr. 601.1-M 90 and UL std N0 60601-1.

7.2.2 Sil.Air 50 TDC compressor

Art. No.	14588
Mains voltage	230 V
Mains frequency	50 / 60 Hz
Mains fuse	T 2,4 A L (SB)
Compressed air output	8 bar
Power consumption	340 VA
Ambient temperature during operation	10-30°C
Ambient temperature during storage	-10 - (+)40°C
Ambient air pressure	500 hPa - 1060 hPa
Air humidity	5 - 95% , non condensing
Weight	22 kg
Housing dimensions (wxhxd)	380 x 330 x 290 mm

7.3 Conformity with standards

EN 60601-1	
- Type of protection against electric shocks:	class 1
- Degree of protection against electric shocks:	type B applied part

7.4 Conformity with directives

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

CE 1275

8 Warranty and Service

8.1 Warranty

During the two years' warranty period from the date of delivery of the product to the end customer, defects will be remedied at no charge to the customer upon the customer furnishing adequate proof that the defect is due to defects in material or workmanship. Shipping costs and risk of loss of returned products shall be borne by the customer.

Please complete the attached warranty card and return it as soon as possible to the address below:

STORZ MEDICAL AG
Unterseestrasse 47
8280 Kreuzlingen
Switzerland

ATTENTION

Any unauthorized opening, repair or modification of the system will relieve the manufacturer of its liability and responsibility for safe system operation and will automatically void the warranty even before the end of the warranty period.



8.2 Service

Should you have any further questions or require additional information, please feel free to call your dealer.